

Biotin anti-human TSLP

Goat Polyclonal Ig
Poly5159



BioLegend
From Revvity

50 µg
Artikel-Nr.: B515903 | BioLegend | Hersteller-Nr.: 515903

300,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: ELISA Detection
Clone: Poly5159
Format: BIOTIN
Isotype: Goat Polyclonal Ig
Species of Origin: Goat
Verpackung: 50 µg

Biotin anti-human TSLP

Clone: Poly5159
Isotype: Goat Polyclonal Ig
Other_Names: Thymic stromal lymphopoietin (TSLP)
Size: 50 µg

Reactivity: Human
Format: BIOTIN
App_Abbrev: ELISA Detection
Species of Origin: Goat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=8150>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>