

APC anti-human CD70

Mouse IgG1, kappa
113-16



BioLegend
From Revvity

100 tests

Artikel-Nr.: B355110 | BioLegend | Hersteller-Nr.: 355110

290,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC
Clone: 113-16
Format: APC
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

APC anti-human CD70

Clone: 113-16

Isotype: Mouse IgG1, kappa

Other_Names: CD27 ligand, CD27L, Ki-24 antigen, Tumor necrosis factor ligand superfamily, member 7 (TNFSF7)

Size: 100 tests

Reactivity: Human

Format: APC

App_Abbrev: FC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=8852>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>