

LEGEND MAX Human CCL17 (TARC) ELISA Kit



BioLegend
From Revvity

1 Pre-coated Plate

Artikel-Nr.: B441107 | BioLegend | Hersteller-Nr.: 441107

370,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: ELISA
Format: LEGEND MAX
Species of Origin: N/A
Verpackung: 1 Pre-coated Plate

LEGEND MAX Human CCL17 (TARC) ELISA Kit

Clone:
Isotype:
Other_Names: Human CCL17 (TARC) Precoated ELISA Kit, Thymus and activation-regulated chemokine (TARC), Small-inducible cytokine 17 (SCYA17), ABCD-2
Size: 1 Pre-coated Plate

Reactivity: Human
Format: LEGEND MAX
App_Abbrev: ELISA
Species of Origin: N/A

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=8865>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>