

PerCP/Cyanine5.5 anti-human RANTES (CCL5)

Mouse IgG2b, kappa
VL1



BioLegend
From Revvity

25 tests

Artikel-Nr.: B515507 | [BioLegend](#) | Hersteller-Nr.: 515507

141,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: VL1
Format: PERCP_CYANINE5.5
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 25 tests

PerCP/Cyanine5.5 anti-human RANTES (CCL5)

Clone: VL1

Isotype: Mouse IgG2b, kappa

Other_Names: CCL5, C-C motif chemokine 5, Small-inducible cytokine A5, T-cell-specific protein RANTES, SIS-delta, T cell-specific protein P228, TCP228, Eosinophil-chemotactic cytokine, EoCP

Size: 25 tests

Reactivity: Human

Format: PERCP_CYANINE5.5

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=8981>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>