

Biotin anti-human CCL17 (TARC)

Goat Polyclonal Ig
Poly5230



BioLegend
From Revvity

50 µg
Artikel-Nr.: B523003 | BioLegend | Hersteller-Nr.: 523003

300,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: ELISA Detection
Clone: Poly5230
Format: BIOTIN
Isotype: Goat Polyclonal Ig
Species of Origin: Goat
Verpackung: 50 µg

Biotin anti-human CCL17 (TARC)

Clone: Poly5230

Isotype: Goat Polyclonal Ig

Other_Names: Thymus and activation-regulated chemokine (TARC), Small-inducible cytokine 17 (SCYA17), ABCD-2

Size: 50 µg

Reactivity: Human

Format: BIOTIN

App_Abbrev: ELISA Detection

Species of Origin: Goat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=9169>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>