

## PE anti-human RANTES (CCL5)

Mouse IgG2b, kappa  
VL1



**BioLegend**  
From Revvity

25 tests

Artikel-Nr.: B515503 | BioLegend | Hersteller-Nr.: 515503

**104,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** ICFC  
**Clone:** VL1  
**Format:** PE  
**Isotype:** Mouse IgG2b, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 25 tests

PE anti-human RANTES (CCL5)

Clone: VL1

Isotype: Mouse IgG2b, kappa

Other\_Names: CCL5, C-C motif chemokine 5, Small-inducible cytokine A5, T-cell-specific protein RANTES, SIS-delta, T cell-specific protein P228, TCP228, Eosinophil-chemotactic cytokine, EoCP

Size: 25 tests

Reactivity: Human

Format: PE

App\_Abbrev: ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=6112>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>