

Recombinant Human LIF (carrier-free)



BioLegend
From Revvity

25 µg
Artikel-Nr.: B593904 | BioLegend | Hersteller-Nr.: 593904

179,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev:BA
Format: RECOM
Verpackung: 25 µg

Recombinant Human LIF (carrier-free)

Clone:

Isotype:

Other_Names: Leukemia inhibitory factor, HILDA (human interleukin in DA cells), MLPLI (melanoma-derived LPL inhibitor), DIA (differentiation-inducing factor)

Size: 25 µg

Reactivity: Human

Format: RECOM

App_Abbrev: BA

Species of Origin:

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=9309>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>