

## Recombinant Human LIF (carrier-free)



**BioLegend**  
From Revvity

100 µg

Artikel-Nr.: B593906 | BioLegend | Hersteller-Nr.: 593906

**615,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

App\_Abbrev:BA

Format: RECOM

Verpackung: 100 µg

Recombinant Human LIF (carrier-free)

Clone:

Isotype:

Other\_Names: Leukemia inhibitory factor, HILDA (human interleukin in DA cells), MLPLI (melanoma-derived LPL inhibitor), DIA (differentiation-inducing factor)

Size: 100 µg

Reactivity: Human

Format: RECOM

App\_Abbrev: BA

Species of Origin:

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=9309>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>