

## Recombinant Mouse LIF (carrier-free)



**BioLegend**  
From Revvity

10 µg  
Artikel-Nr.: B554002 | BioLegend | Hersteller-Nr.: 554002

**113,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

App\_Abbrev:BA  
Format: RECOM  
Verpackung: 10 µg

Recombinant Mouse LIF (carrier-free)

Clone:  
Isotype:  
Other\_Names: MLPLI (melanoma-derived LPL inhibitor), DIA (differentiation-inducing factor), LIF (leukemia inhibitory factor)  
Size: 10 µg

Reactivity: Mouse  
Format: RECOM  
App\_Abbrev: BA  
Species of Origin:

Weitere Informationen zum Produkt. Link zum Hersteller:  
<https://www.biolegend.com/Default.aspx?ID=6664&productid=10032>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>