

Recombinant Mouse LIF (carrier-free)



BioLegend
From Revvity

25 µg
Artikel-Nr.: B554004 | BioLegend | Hersteller-Nr.: 554004

179,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev:BA
Format: RECOM
Verpackung: 25 µg

Recombinant Mouse LIF (carrier-free)

Clone:
Isotype:
Other_Names: MLPLI (melanoma-derived LPL inhibitor), DIA (differentiation-inducing factor), LIF (leukemia inhibitory factor)
Size: 25 µg

Reactivity: Mouse
Format: RECOM
App_Abbrev: BA
Species of Origin:

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=10032>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>