

Recombinant Mouse LIF (carrier-free)



BioLegend
From Revvity

500 µg
Artikel-Nr.: B554008 | **BioLegend** | Hersteller-Nr.: 554008

1.496,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev:FC, ICFC

Format: RECOM

Verpackung: 500 µg

Recombinant Mouse LIF (carrier-free)

Clone:

Isotype:

Other_Names: MLPLI (melanoma-derived LPL inhibitor), DIA (differentiation-inducing factor), LIF (leukemia inhibitory factor)

Size: 500 µg

Reactivity: Mouse

Format: RECOM

App_Abbrev: FC, ICFC

Species of Origin:

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=10032>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>