

Purified anti-human CCL5 (RANTES)

Mouse IgG2b, kappa
J047C5



BioLegend
From Revvity

500 µg

Artikel-Nr.: B526402 | BioLegend | Hersteller-Nr.: 526402

278,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: ELISA Capture
Clone: J047C5
Format: PURE
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 500 µg

Purified anti-human CCL5 (RANTES)

Clone: J047C5

Isotype: Mouse IgG2b, kappa

Other_Names: SIS-delta, T-cell specific protein p288, beta-chemokine RANTES, regulated upon activation, normally T-expressed, and presumably secreted, small inducible cytokine A5

Size: 500 µg

Reactivity: Human

Format: PURE

App_Abbrev: ELISA Capture

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=10768>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>