

Biotin anti-human CCL5 (RANTES)

Goat Polyclonal Ig
Poly5197



BioLegend
From Revvity

50 µg
Artikel-Nr.: B519703 | BioLegend | Hersteller-Nr.: 519703

300,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: ELISA Detection
Clone: Poly5197
Format: BIOTIN
Isotype: Goat Polyclonal Ig
Species of Origin: Goat
Verpackung: 50 µg

Biotin anti-human CCL5 (RANTES)

Clone: Poly5197

Isotype: Goat Polyclonal Ig

Other_Names: SIS-delta, T-cell specific protein p288, beta-chemokine RANTES, regulated upon activation, normally T-expressed, and presumably secreted, small inducible cytokine A5

Size: 50 µg

Reactivity: Human

Format: BIOTIN

App_Abbrev: ELISA Detection

Species of Origin: Goat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=10633>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>