

## Purified anti-V5 Epitope Tag

Rabbit Polyclonal IgG  
Poly29038



BioLegend  
From Revvity

100 µL

Artikel-Nr.: B903801 | BioLegend | Hersteller-Nr.: 903801

**265,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** ELISA, ICC, IHC, Purification, WB  
**Clone:** Poly29038  
**Format:** PURE  
**Isotype:** Rabbit Polyclonal IgG  
**Species of Origin:** Rabbit  
**Verpackung:** 100 µl

Purified anti-V5 Epitope Tag

Clone: Poly29038

Isotype: Rabbit Polyclonal IgG

Other\_Names: GKPIPPLLGLDST epitope tag, GKPIPPLLGLDST tag, Protein Rev, Regulator of expression of viral proteins, Rev, V5 epitope tag

Size: 100 µL

Reactivity:

Format: PURE

App\_Abbrev: WB, ICC, IHC, ELISA, Purification

Species of Origin: Rabbit

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=11593>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>