

Purified anti-human IgG1 Fc

Rat IgG2a, kappa
M1308A10



BioLegend
From Revvity

100 µg

Artikel-Nr.: B410601 | BioLegend | Hersteller-Nr.: 410601

90,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: WB
Clone: M1308A10
Format: PURE
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 100 µg

Purified anti-human IgG1 Fc

Clone: M1308A10
Isotype: Rat IgG2a, kappa
Other_Names: Immunoglobulin Heavy Constant Gamma 1
Size: 100 µg

Reactivity: Human
Format: PURE
App_Abbrev: WB
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=10766>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>