

Purified anti-E1 Ubiquitin Activating Enzyme

Mouse IgG1, kappa
2G2.3.5



BioLegend
From Revvity

100 µg

Artikel-Nr.: B932002 | BioLegend | Hersteller-Nr.: 932002

249,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ICC, WB
Clone: 2G2.3.5
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-E1 Ubiquitin Activating Enzyme

Clone: 2G2.3.5

Isotype: Mouse IgG1, kappa

Other_Names: Ubiquitin-Activating Enzyme E1, Ubiquitin-Like Modifier Activating Enzyme 1, UBE1, A1S9T, GXP1, SMAX2, 2G.2.3.5, 2G2.3-5

Size: 100 µg

Reactivity: Human, Mouse

Format: PURE

App_Abbrev: WB, ICC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=11850>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>