

Alexa Fluor® 647 anti-mouse CD20

Rat IgG2b, kappa
SA275A11



BioLegend
From Revvity

25 µg
Artikel-Nr.: B150403 | BioLegend | Hersteller-Nr.: 150403

128,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, IHC-F
Clone: SA275A11
Format: A647
Isotype: Rat IgG2b, kappa
Species of Origin: Rat
Verpackung: 25 µg

Alexa Fluor® 647 anti-mouse CD20

Clone: SA275A11
Isotype: Rat IgG2b, kappa
Other_Names: Ms4a1, Ly-44
Size: 25 µg

Reactivity: Mouse
Format: A647
App_Abbrev: FC, IHC-F
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=12238>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>