

Purified anti-Arginase 1

Mouse IgG1, kappa
O94E6/ARG1



BioLegend
From Revvity

25 µg
Artikel-Nr.: B678801 | BioLegend | Hersteller-Nr.: 678801

81,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: IHC-P, KO/KD-WB, WB
Clone: O94E6/ARG1
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-Arginase 1

Clone: O94E6/ARG1
Isotype: Mouse IgG1, kappa
Other_Names: Arginase 1, liver type Arginase, Type I Arginase, ARG1
Size: 25 µg

Reactivity: Human, Mouse
Format: PURE
App_Abbrev: WB, KO/KD-WB, IHC-P
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=12106>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>