

Purified anti-SNAP-25

Mouse IgG1, kappa
SMI 81



BioLegend
From Revvity

100 µg

Artikel-Nr.: B836304 | BioLegend | Hersteller-Nr.: 836304

221,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: EM, FC, ICC, IHC-P, IP, WB
Clone: SMI 81
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-SNAP-25

Clone: SMI 81

Isotype: Mouse IgG1, kappa

Other_Names: Synaptosomal-associated protein 25, SUP, Synaptosomal-associated 25 kD Protein

Size: 100 µg

Reactivity: Human, Mouse, Rat

Format: PURE

App_Abbrev: IHC-P, WB, EM, FC, ICC, IP

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=12727>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>