

Purified anti-human LIF

Goat Polyclonal IgG
Poly5233



BioLegend
From Revvity

100 µg
Artikel-Nr.: B523301 | **BioLegend** | Hersteller-Nr.: 523301

280,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: Poly5233
Format: PURE
Isotype: Goat Polyclonal IgG
Species of Origin: Goat
Verpackung: 100 µg

Purified anti-human LIF

Clone: Poly5233

Isotype: Goat Polyclonal IgG

Other_Names: Leukemia Inhibitory Factor, HILDA (human interleukin in DA cells), MLPLI (melanoma-derived LPL inhibitor), DIA (differentiation-inducing factor)

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Goat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=12804>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>