

Biotin anti-human LIF

Mouse IgG1, kappa
M017C3



BioLegend
From Revvity

50 µg
Artikel-Nr.: B530301 | BioLegend | Hersteller-Nr.: 530301

141,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ELISA Detection
Clone: M017C3
Format: BIOTIN
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 50 µg

Biotin anti-human LIF

Clone: M017C3

Isotype: Mouse IgG1, kappa

Other_Names: Leukemia Inhibitory Factor, HILDA (human interleukin in DA cells), MLPLI (melanoma-derived LPL inhibitor), DIA (differentiation-inducing factor)

Size: 50 µg

Reactivity: Human

Format: BIOTIN

App_Abbrev: ELISA Detection

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=12807>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>