

## Purified anti-human Arginase I

Mouse IgG2b, kappa  
14D2C43



**BioLegend**  
From Revvity

25 µg  
Artikel-Nr.: B369701 | BioLegend | Hersteller-Nr.: 369701

**81,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** ICFC  
**Clone:** 14D2C43  
**Format:** PURE  
**Isotype:** Mouse IgG2b, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 25 µg

Purified anti-human Arginase I

Clone: 14D2C43  
Isotype: Mouse IgG2b, kappa  
Other\_Names: ARG1, Liver Arginase, Type 1 Arginase, Arginase-1, Liver-Type Arginase  
Size: 25 µg

Reactivity: Human  
Format: PURE  
App\_Abbrev: ICFC  
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:  
<https://www.biolegend.com/Default.aspx?ID=6664&productid=12911>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>