

LEGEND MAX Mouse LIF ELISA Kit



BioLegend
From Revvity

1 Pre-coated Plate
Artikel-Nr.: B445107 | BioLegend | Hersteller-Nr.: 445107

398,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: ELISA
Format: LEGEND MAX
Species of Origin:N/A
Verpackung: 1 Pre-coated Plate

LEGEND MAX Mouse LIF ELISA Kit

Clone:
Isotype:
Other_Names: Leukemia inhibitory factor, CDF, DIA, HILDA, MLPLI
Size: 1 Pre-coated Plate

Reactivity: Mouse
Format: LEGEND MAX
App_Abbrev: ELISA
Species of Origin: N/A

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=13433>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>