

Hot Start Taq DNA Polymerase, 5 U/ μ l



500 Units

Artikel-Nr.: 350200102 | biotechrabbit | Hersteller-Nr.: BR0200102

93,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

High-Fidelity:	nein
Hot-Start:	ja
Master Mix / Enzym:	Enzym
Multiplex:	nein
PCR / RTPCR:	PCR
PCR-Optimierung:	Standard
Verpackung:	500 Units