

MAGic Clamp COMBO

for 125 ml, 250 ml, 500 ml flasks



55H1000-MR-CMB Adjustable flask clamp, accepts 125ml, 250ml and 500ml Erlenmeyer flasks.
(For a video demonstration, scan the QR code above.)

Benchmark
Scientific

1 Stück
Artikel-Nr.: 55H1000-MR-CMB | Benchmark | Hersteller-Nr.: H1000-MR-CMB

58,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

Magic Clamp: ja

Verpackung: 1 Stück