

## Goat anti-rabbit second. antibody AzureSpectra 650

AzureSpectra 650 conjugated secondary antibody  
for fluorescent Westerns



250 µl

Artikel-Nr.: 512165 | Azure

**125,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**Anregung:** 653 nm  
**Anti-:** Rabbit  
**Applikation:** Fluoreszenz  
**Emission:** 672 nm  
**Konjugat:** Fluorophor  
**Origin:** Goat  
**Verpackung:** 250 µl