

## Purified anti-HMB-45

Mouse IgG1, kappa  
HMB-45



**BioLegend**  
From Revvity

100 µg

Artikel-Nr.: B911504 | BioLegend | Hersteller-Nr.: 911504

**249,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** IHC-F, IHC-P  
**Clone:** HMB-45  
**Format:** PURE  
**Isotype:** Mouse IgG1, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 100 µg

Purified anti-HMB-45

Clone: HMB-45

Isotype: Mouse IgG1, kappa

Other\_Names: Melanoma Marker antibody, Human Melanoma Black, Melanocyte protein PMEL, PMEL

Size: 100 µg

Reactivity: Human

Format: PURE

App\_Abbrev: IHC-P, IHC-F

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=14906>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>