

AzureSpectra Labeling Kit 700

Kit for covalently labeling antibodies with a fluorescent dye



1 Stück
Artikel-Nr.: 512188 | Azure

517,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

Verpackung: 1 Stück

AzureSpectra Labeling Kit 700 - Kit for covalently labeling antibodies with a fluorescent dye – includes labeling buffer, dye solution, quenching solution, neutralization buffer and two G-25 columns – each kit includes sufficient reagents to perform two labelings of 50 micrograms of antibody