

## AzureSpectra Labeling Kit 800

Kit for covalently labeling antibodies with a fluorescent dye



1 Stück  
Artikel-Nr.: 512189 | Azure

**517,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**Verpackung:** 1 Stück

AzureSpectra Labeling Kit 800 - Kit for covalently labeling antibodies with a fluorescent dye – includes labeling buffer, dye solution, quenching solution, neutralization buffer and two G-25 columns – each kit includes sufficient reagents to perform two labelings of 50 micrograms of antibody