

## Donkey anti-goat secondary AB AzureSpectra 550

AzureSpectra 550 conjugated secondary antibody  
for fluorescent Westerns



250 µl

Artikel-Nr.: 512164 | Azure

**125,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**Anregung:** 551 nm  
**Anti-:** Goat  
**Applikation:** Fluoreszenz  
**Emission:** 565 nm  
**Konjugat:** Fluorophor  
**Origin:** Donkey  
**Verpackung:** 250 µl