

PerCP/Cyanine5.5 anti-mouse/human CD324 (E-Cad...

Rat IgG1, kappa
DECMA-1



BioLegend
From Revvity

100 µg

Artikel-Nr.: B147318 | **BioLegend** | Hersteller-Nr.: 147318

306,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC
Clone: DECMA-1
Format: PERCP_CYANINE5.5
Isotype: Rat IgG1, kappa
Species of Origin: Rat
Verpackung: 100 µg

PerCP/Cyanine5.5 anti-mouse/human CD324 (E-Cadherin)

Clone: DECMA-1
Isotype: Rat IgG1, kappa
Other_Names: E-Cadherin, Cadherin-1, CDH1, and UVO
Size: 100 µg

Reactivity: Mouse, Human
Format: PERCP_CYANINE5.5
App_Abbrev: FC
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=16389>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>