

Alexa Fluor® 488 anti-human CD27 Recombinant

Mouse IgG1, kappa
QA17A18



BioLegend
From Revvity

25 tests

Artikel-Nr.: B393203 | BioLegend | Hersteller-Nr.: 393203

95,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC
Clone: QA17A18
Format: A488
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 tests

Alexa Fluor® 488 anti-human CD27 Recombinant

Clone: QA17A18
Isotype: Mouse IgG1, kappa
Other_Names: S152, T14, TNFRSF7
Size: 25 tests

Reactivity: Human
Format: A488
App_Abbrev: FC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=16257>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>