

Biotin anti-mouse LIF

Goat Polyclonal IgG
Poly5232



BioLegend
From Revvity

50 µg
Artikel-Nr.: B523204 | BioLegend | Hersteller-Nr.: 523204

300,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: ELISA Detection
Clone: Poly5232
Format: BIOTIN
Isotype: Goat Polyclonal IgG
Species of Origin: Goat
Verpackung: 50 µg

Biotin anti-mouse LIF

Clone: Poly5232
Isotype: Goat Polyclonal IgG
Other_Names: Leukemia Inhibitory factor, HILDA, CDF, DIA, MLPLI
Size: 50 µg

Reactivity: Mouse
Format: BIOTIN
App_Abbrev: ELISA Detection
Species of Origin: Goat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=17973>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>