

Purified anti-Aldose Reductase (AKR1B1)

Mouse IgG1, kappa
4G9-B6-F6



BioLegend
From Revvity

25 µg
Artikel-Nr.: B870301 | BioLegend | Hersteller-Nr.: 870301

85,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: IHC-P, WB
Clone: 4G9-B6-F6
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-Aldose Reductase (AKR1B1)

Clone: 4G9-B6-F6

Isotype: Mouse IgG1, kappa

Other_Names: Aldo-keto reductase family 1 member B1, Aldehyde reductase, Aldose reductase, AR

Size: 25 µg

Reactivity: Human, Mouse

Format: PURE

App_Abbrev: WB, IHC-P

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=18678>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>