

## Purified anti-Aldose Reductase (AKR1B1)

Mouse IgG1, kappa  
4G9-B6-F6



**BioLegend**  
From Revvity

100 µg  
Artikel-Nr.: B870302 | BioLegend | Hersteller-Nr.: 870302

**221,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** IHC-P, WB  
**Clone:** 4G9-B6-F6  
**Format:** PURE  
**Isotype:** Mouse IgG1, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 100 µg

Purified anti-Aldose Reductase (AKR1B1)

Clone: 4G9-B6-F6  
Isotype: Mouse IgG1, kappa  
Other\_Names: Aldo-keto reductase family 1 member B1, Aldehyde reductase, Aldose reductase, AR  
Size: 100 µg

Reactivity: Human, Mouse  
Format: PURE  
App\_Abbrev: WB, IHC-P  
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:  
<https://www.biolegend.com/Default.aspx?ID=6664&productid=18678>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>