

Purified anti-CDC20

Rat IgG2a, kappa
W18156A



BioLegend
From Revvity

25 µg
Artikel-Nr.: B934701 | **BioLegend** | Hersteller-Nr.: 934701

81,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: W18156A
Format: PURE
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 25 µg

Purified anti-CDC20

Clone: W18156A

Isotype: Rat IgG2a, kappa

Other_Names: Cell Division Cycle 20, Cell Division Cycle Protein 20 Homolog, BA276H19.3, CDC20A, BA276H19.3, p55CDC

Size: 25 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=18652>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>