

Filterspitzen 20 µl BevelPoint

steril
Nachfolgeartikel für Art.-Nr. 701041



labcon

10 Racks à 96 Stück
Artikel-Nr.: LC2555 | Labcon | Hersteller-Nr.: 1055-965-018-9

77,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

Farbe: farblos
Filter: ja
Länge: normal
Low Binding: nein
Maximalvolumen: 20 µl
Pipettenspitzen Linie: Labcon
Sterilität: steril
Verpackung: 10 Racks à 96 Stück

Filterspitze passend für:

- Brand Transferpette 2-20 µl, 8/12 Kanal 2-20 µl, 2.5-25 µl
- Eppendorf Reference 2 2-20 µl
- Eppendorf Research Pro 2-20 µl
- Eppendorf Xplorer (plus) 2-20 µl
- Gilson Neo 2-20 µl
- Gilson Pipetman Classic G, L 2-20 µl
- Gilson Pipetman M 2-20 µl, 8/12 Kanal 1-20 µl
- Gilson Pipetman Ultra 2-20 µl, 8/12 Kanal 1-20 µl
- Nichiryo 20 µl
- Rainin Classic 2-20 µl
- Rainin EDP1 2-20 µl

- Rainin EDP3 Plus 2-20 µl
- Rainin Pipet Lite SL 2-20 µl
- Sartorius mLine 2-20 µl
- Sartorius Proline plus, Tacta 2-20 µl
- Socorex Acura 825, Manual XS 2-20 µl
- Thermo Digital, F1, F2 2-20 µl
- Thermo Focus 2-20 µl