

## TotalSeq-B0902 anti-human CD328 (Siglec-7)

Mouse IgG1, kappa  
6-434



**BioLegend**  
From Revvity

10 µg  
Artikel-Nr.: B339219 | BioLegend | Hersteller-Nr.: 339219

**296,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** PG  
**Clone:** 6-434  
**Format:** TOTALSEQ\_B  
**Isotype:** Mouse IgG1, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 10 µg

TotalSeq-B0902 anti-human CD328 (Siglec-7)

Clone: 6-434  
Isotype: Mouse IgG1, kappa  
Other\_Names: Siglec7 (Sialic-binding Ig like lectin 7), p75, AIRM1  
Size: 10 µg

Reactivity: Human  
Format: TOTALSEQ\_B  
App\_Abbrev: PG  
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:  
<https://www.biolegend.com/Default.aspx?ID=6664&productid=20067>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>