

LEAF Purified anti-human VEGF-165

Goat Polyclonal IgG
Poly5225



BioLegend
From Revvity

100 µg

Artikel-Nr.: B522505 | BioLegend | Hersteller-Nr.: 522505

362,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: Block, WB
Clone: Poly5225
Format: LEAF
Isotype: Goat Polyclonal IgG
Species of Origin: Goat
Verpackung: 100 µg

LEAF Purified anti-human VEGF-165

Clone: Poly5225

Isotype: Goat Polyclonal IgG

Other_Names: VEGFA165, MVCD1, VEGF, Vascular permeability factor (VPF)

Size: 100 µg

Reactivity: Human

Format: LEAF

App_Abbrev: Block, WB

Species of Origin: Goat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=21056>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>