

Ultra-LEAF Purified anti-human CXCL9 (MIG)

Mouse IgG1, kappa
A16023A



BioLegend
From Revvity

1 mg

Artikel-Nr.: B943504 | BioLegend | Hersteller-Nr.: 943504

836,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: Neut
Clone: A16023A
Format: ULEAF
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 1 mg

Ultra-LEAF Purified anti-human CXCL9 (MIG)

Clone: A16023A

Isotype: Mouse IgG1, kappa

Other_Names: Monokine induced by interferon gamma (MIG), Small inducible cytokine, subfamily B, member 9 (SCYB9)

Size: 1 mg

Reactivity: Human

Format: ULEAF

App_Abbrev: Neut

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=20949>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>