

## Ultra-LEAF Purified anti-human CXCL11 (ITAC)

Mouse IgG2a, kappa  
A15129H



BioLegend  
From Revvity

100 µg

Artikel-Nr.: B945403 | BioLegend | Hersteller-Nr.: 945403

**287,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** Neut  
**Clone:** A15129H  
**Format:** ULEAF  
**Isotype:** Mouse IgG2a, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 100 µg

Ultra-LEAF Purified anti-human CXCL11 (ITAC)

Clone: A15129H

Isotype: Mouse IgG2a, kappa

Other\_Names: Interferon-inducible T-cell alpha chemoattractant (ITAC), Small inducible cytokine, subfamily B, member 9B (SCYB9B), Small inducible cytokine, subfamily B, member 11 (SCYB11)

Size: 100 µg

Reactivity: Human

Format: ULEAF

App\_Abbrev: Neut

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=20893>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>