

Ultra-LEAF Purified anti-human CXCL11 (ITAC)

Mouse IgG2a, kappa
A15129H



BioLegend
From Revvity

1 mg
Artikel-Nr.: B945404 | BioLegend | Hersteller-Nr.: 945404

964,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: Neut
Clone: A15129H
Format: ULEAF
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 1 mg

Ultra-LEAF Purified anti-human CXCL11 (ITAC)

Clone: A15129H

Isotype: Mouse IgG2a, kappa

Other_Names: Interferon-inducible T-cell alpha chemoattractant (ITAC), Small inducible cytokine, subfamily B, member 9B (SCYB9B), Small inducible cytokine, subfamily B, member 11 (SCYB11)

Size: 1 mg

Reactivity: Human

Format: ULEAF

App_Abbrev: Neut

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=20893>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>