

PE/Cyanine5 anti-human CD326 (EpCAM)

Mouse IgG2b, kappa
9C4



BioLegend
From Revvity

25 tests

Artikel-Nr.: B324251 | BioLegend | Hersteller-Nr.: 324251

147,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC
Clone: 9C4
Format: PE_CYANINE5
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 25 tests

PE/Cyanine5 anti-human CD326 (EpCAM)

Clone: 9C4

Isotype: Mouse IgG2b, kappa

Other_Names: Ep-CAM, tumor associated calcium signal transducer 1, epithelial cell surface antigen, epithelial glycoprotein 2, EGP2, adenocarcinoma associated antigen, TROP1

Size: 25 tests

Reactivity: Human

Format: PE_CYANINE5

App_Abbrev: FC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=21147>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>