

Ultra-LEAF Purified anti-human VEGF-165

Rat IgG2b, kappa
A20088B



BioLegend
From Revvity

100 µg

Artikel-Nr.: B668053 | BioLegend | Hersteller-Nr.: 668053

331,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: Block
Clone: A20088B
Format: ULEAF
Isotype: Rat IgG2b, kappa
Species of Origin: Rat
Verpackung: 100 µg

Ultra-LEAF Purified anti-human VEGF-165

Clone: A20088B
Isotype: Rat IgG2b, kappa
Other_Names: VEGFA165, MVCD1, VEGF, Vascular permeability factor (VPF)
Size: 100 µg

Reactivity: Human
Format: ULEAF
App_Abbrev: Block
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=21885>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>