

LEAF Purified anti-mouse SCF

Goat Polyclonal IgG
Poly5417



BioLegend
From Revvity

100 µg

Artikel-Nr.: B541703 | BioLegend | Hersteller-Nr.: 541703

288,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: Block, WB
Clone: Poly5417
Format: LEAF
Isotype: Goat Polyclonal IgG
Species of Origin: Goat
Verpackung: 100 µg

LEAF Purified anti-mouse SCF

Clone: Poly5417

Isotype: Goat Polyclonal IgG

Other_Names: Kit Ligand, KITLG, Mast cell growth factor, MGF, Steel factor, SLF, Stem Cell Factor

Size: 100 µg

Reactivity: Mouse

Format: LEAF

App_Abbrev: Block, WB

Species of Origin: Goat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=22807>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>