

Purified anti-Peroxiredoxin 1

Mouse IgG1, kappa
6A6-4H3-7A3



BioLegend
From Revvity

25 µg
Artikel-Nr.: B611151 | BioLegend | Hersteller-Nr.: 611151

96,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ICC, WB
Clone: 6A6-4H3-7A3
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-Peroxiredoxin 1

Clone: 6A6-4H3-7A3

Isotype: Mouse IgG1, kappa

Other_Names: PRDX1, Natural killer cell-enhancing factor A, NKEF-A, Proliferation-associated gene protein, PAG, Thioredoxin peroxidase 2, Thioredoxin-dependent peroxide reductase 2, Thioredoxin-dependent peroxiredoxin 1, PAGA, PAGB, TDPX2

Size: 25 µg

Reactivity: Human, Mouse, Rat

Format: PURE

App_Abbrev: WB, ICC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=22833>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>