

Caretta 96-Well PCR Platten, farblos, barcode

PC Full Skirt, Low Profile, Barcode, Eckschnitt H1
DNA-, DNase-, RNasefrei



Packung à 100 Stück
Artikel-Nr.: 711810 | Biozym

703,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Cyclertyp: Analytic Jena / Biometra-FlexCycler2, Analytic Jena / Biometra-T1, Analytic Jena / Biometra-TAdvanced 96, Analytic Jena / Biometra-TAdvanced 96G, Analytic Jena / Biometra-TAdvanced 96S, Analytic Jena / Biometra-TAdvanced 96SG, Analytic Jena / Biometra-TGradient, Analytic Jena / Biometra-TOne, Analytic Jena / Biometra-TProfessional Serie, Analytic Jena / Biometra-qTOWERiris (touch), Bio-Rad / MJ Research-C1000 (Touch), Bio-Rad / MJ Research-CFX96 Touch Deep Well, Bio-Rad / MJ Research-DNA Engine Tetrad (2), Bio-Rad / MJ Research-DNA Engine, Bio-Rad / MJ Research-DNA Engine Dyad, Bio-Rad / MJ Research-DNA Engine Dyad Disciple, Bio-Rad / MJ Research-PC Tempo 96 / Deep Well, Bio-Rad / MJ Research-PTC 200, Bio-Rad / MJ Research-S1000, Bioer-GeneTouch 96, Bioer-GeneTouch Plus 96, Bioer-LifeTouch, Bioer-Lifeeco, Eppendorf-Mastercycler ep (gradient), Eppendorf-Mastercycler ep pro (S), Eppendorf-Mastercycler ep realplex, Eppendorf-Mastercycler nexus (gradient), Eppendorf-Mastercycler nexus GSX1, Eppendorf-Mastercycler nexus SX1, VWR-Collection UNO96

Eckschnitt: H1

Eigenschaften: mit Barcode, rigid frame

Farbe: Rahmen: farblos, farblos

Format: full skirted

Kategorie: 96-Well Platte

Produktlinie: Caretta

Profil: Low Profile

Verpackung: Packung à 100 Stück

Volumen: 0.1 ml