

LEAF Purified anti-mouse LIF

Goat Polyclonal IgG
Poly5232



BioLegend
From Revvity

100 µg

Artikel-Nr.: B523205 | **BioLegend** | Hersteller-Nr.: 523205

320,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: Poly5232
Format: LEAF
Isotype: Goat Polyclonal IgG
Species of Origin: Goat
Verpackung: 100 µg

LEAF Purified anti-mouse LIF

Clone: Poly5232
Isotype: Goat Polyclonal IgG
Other_Names: Leukemia Inhibitory factor, HILDA, CDF, DIA, MLPLI
Size: 100 µg

Reactivity: Mouse
Format: LEAF
App_Abbrev: FC, ICFC
Species of Origin: Goat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=23187>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>