

## Purified anti-human HGF

Mouse IgG2a, kappa  
W19372J



**BioLegend**  
From Revvity

100 µg

Artikel-Nr.: B613052 | BioLegend | Hersteller-Nr.: 613052

**188,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** Direct ELISA  
**Clone:** W19372J  
**Format:** PURE  
**Isotype:** Mouse IgG2a, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 100 µg

Purified anti-human HGF

Clone: W19372J

Isotype: Mouse IgG2a, kappa

Other\_Names: Scatter Factor (SF), Hepatopietin (HPTA), Lung fibroblast-derived mitogen

Size: 100 µg

Reactivity: Human

Format: PURE

App\_Abbrev: Direct ELISA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=23115>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>