

Spark Red 718 anti-human CD69

Mouse IgG1, kappa
FN50



BioLegend
From Revvity

100 tests

Artikel-Nr.: B310968 | **BioLegend** | Hersteller-Nr.: 310968

196,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: FN50
Format: SPKR_718
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

Spark Red 718 anti-human CD69

Clone: FN50

Isotype: Mouse IgG1, kappa

Other_Names: Very Early Activation Antigen (VEA), Activation inducer molecule (AIM)

Size: 100 tests

Reactivity: Human

Format: SPKR_718

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=23489>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>